

Know Cancer's Next Move

- Does my patient need adjuvant chemotherapy?
- Is the treatment working?
- Is the cancer recurring?

Discover unique insights with molecular residual disease (MRD) detection, informed by ctDNA

AUTHORIZED DISTRIBUTOR



Signatera™
Residual disease test (MRD)

The LifeLabs logo features the word "LifeLabs" in a large, blue, sans-serif font, with a red dot above the "i". Below it, the word "GENETICS" is written in a smaller, blue, sans-serif font, with a red dot above the "I".


LifeLabs
GENETICS®

When to use Signatera™ for patients with CRC

In the adjuvant setting

Use after surgery to evaluate the need for adjuvant chemotherapy and potentially avoid unnecessary treatment

Signatera™ MRD status outperforms known clinicopathologic risk factors in predicting relapse¹⁻⁴



- > 97% of patients with a positive Signatera™ result relapsed without additional treatment¹
- > Serial testing improved sensitivity and negative predictive value of results


In the surveillance setting

Use along with CEA testing and other surveillance modalities to detect recurrence earlier, to enable surgical resection or other early intervention

Signatera™ was shown to detect relapse more accurately than CEA with clinically meaningful lead times over CT scans¹

- > Get clarity when evaluating patients with indeterminate CEA levels or CT scans
- > Signatera™ facilitates shared decision-making and confident treatment planning

Accuracy in detecting relapse



CEA = carcinoembryonic antigen; CT = computed tomography; ctDNA = circulating-tumour DNA

Average lead time
of ctDNA detection
before CT scan


8.7
MONTHS

Maximum lead of 16.5 months


A large-scale prospective, MRD-guided study evaluated the clinical utility of ctDNA analysis in colorectal cancer (CRC)

MRD-positive CRC patients at 4 weeks post-op benefited significantly from chemotherapy while MRD-negative patients at 4 weeks post-op did not demonstrate any significant trend in treatment benefit⁵

MRD-positive patients at 4 weeks post-op



MRD-negative patients at 4 weeks post-op



Track ctDNA dynamics to enable longitudinal monitoring

- Signatera™ reports presence/absence of ctDNA and ctDNA quantity in terms of MTM/mL for longitudinal assessment

Final Results Summary

Signatera™ Negative




MTM/mL

Not Detected

Mean tumour molecules per mL is calculated based on the mean of ctDNA molecules detected per mL of the patient's plasma

Historical Results



MTM = mean tumour molecules

AUTHORIZED DISTRIBUTOR



Signatera™
Residual disease test (MRD)



Just like no two tumours are alike, Signatera™ is personalized for each patient



Tumour-informed MRD assay for individualized care

- Customized for each patient's unique tumour signature using Whole Exome Sequencing (WES) to target the top clonal mutations



High sensitivity and specificity for accurate MRD assessment

- By only tracking tumour-specific variants, sensitivity is optimized with a LOD down to 0.01% VAF⁶
- Filters out germline and CHIP mutations to reduce background noise and to minimize false positives




Reliable longitudinal monitoring for confident decision-making


- Tracks ctDNA dynamics by MTM/mL to enable longitudinal monitoring with a simple blood draw
- Follows clonal mutations that should persist as the tumour evolves

LOD = limit of detection; CHIP = clonal hematopoiesis of indeterminate potential; VAF = Variant allele frequency; WES = whole exome sequencing


Signatera™ is validated across multiple tumour types^{1,6-8}




88%
sensitivity to relapse
Average lead time 8.7 mos



89%
sensitivity to relapse
Average lead time 9.5 mos



92%
sensitivity to relapse
Average lead time 4.0 mos



100%
sensitivity to relapse
Average lead time 2.8 mos



LifeLabs Genetics
175 Galaxy Blvd., Suite 105,
Toronto, ON M9W 0C9
www.lifelabsgenetics.com/signatera



**To place an order or learn more,
contact our Genetics Team**
Ask.Genetics@lifelabs.com
1-844-GENE HELP (1-844-363-4357)

References

1. Reinert T, Henriksen TV, Christensen E, et al. Analysis of Plasma Cell-Free DNA by Ultradeep Sequencing in Patients With Stages I to III Colorectal Cancer. JAMA Oncol. 2019.
2. Sinicrope FA, Foster NR, Thibodeau SN, et al. DNA Mismatch Repair Status and Colon Cancer Recurrence and Survival in Clinical Trials of 5-Fluorouracil-Based Adjuvant Therapy. J Natl Cancer Inst. 2011;103(11):863-875.
3. Aoyama, Oba K, Honda M, et al. Impact of postoperative complications on the colorectal cancer survival and recurrence: analyses of pooled individual patients' data from three large phase III randomized trials. Cancer Med. 2017;6(7):1573-1580.
4. Yothers G, O'Connell MJ, Lopatin M, et al. Validation of the 12-gene colon cancer recurrence score in NSABP C-07 as a predictor of recurrence in patients with stage II and III colon cancer treated with fluorouracil and leucovorin (FU/LV) and FU/LV plus oxaliplatin. J Clin Oncol. 2013;31(36):4512-4519.
5. Kotani D, et al., Molecular residual disease and efficacy of adjuvant chemotherapy in patients with colorectal cancer, Nature Medicine v29 Issue 1 Jan 2023
6. Coombes RC, Page K, Salari R, et al. Personalized Detection of Circulating Tumour DNA Antedates Breast Cancer Metastatic Recurrence. Clin Cancer Res. 2019;25(14):4255-4263.
7. Abbosh C, Birkbak NJ, Wilson GA, et al. Phylogenetic ctDNA analysis depicts early-stage lung cancer evolution. Nature. 2017;545(7655):446-451.
8. Christensen E, Birkenkamp-Demtroder K, Sethi H, et al. Early Detection of Metastatic Relapse and Monitoring of Therapeutic Efficacy by Ultra-Deep Sequencing of Plasma Cell-Free DNA in Patients With Urothelial Bladder Carcinoma. J Clin Oncol. 2019;37(18):1547-1557.
9. Data on file

AUTHORIZED DISTRIBUTOR



Signatera™
Residual disease test (MRD)

LifeLabs Genetics is an authorized distributor of Natera's Signatera™ in Canada. The content has not been reviewed by nor endorsed by Natera, Inc. LifeLabs® is solely responsible for maintaining content according to Natera partnership guidelines as well as all legal and regulatory requirements in Canada. Natera™ is CAP accredited, ISO 13485 and CLIA certified. © 2025 Natera, Inc. All Rights Reserved. LifeLabs is a registered trademark and LifeLabs Genetics is a trademark of LifeLabs LP. LifeLabs LP © LifeLabs 2025.



LifeLabs Genetics Signatera™ Physician Brochure | VER 1.1 APR 2025